

Validating Thermal Neutron Capture γ -Ray Data using the RPI Gaerttner LINAC Center

Katelyn Cook¹, Yaron Danon¹, Amanda Lewis², Michael Rapp², Adam Daskalakis², Devin Barry², Peter Brain¹, Dominik Fritz¹, Alec Golas¹, Adam Ney¹, Gregory Siemers¹, Sukhjinder Singh¹, and Benjamin Wang¹

¹Gaerttner LINAC Center, Rensselaer Polytechnic Institute, Troy, NY

²Naval Nuclear Laboratory, Schenectady, NY



2024 International Topical Meeting on Nuclear Applications of Accelerators (AccelApp24)
Jefferson Lab in Norfolk, Virginia
March 17-21




Rensselaer

NAVAL NUCLEAR
LABORATORY

linac
The Gaertner LINAC Center




Outline

- **Introduction:** RPI Gaerttner LINAC Center, Neutron Capture Reactions, Project Goals & Motivation
- **Experimental Methods:** RPI Capture γ -Ray Multiplicity Detector
- **Simulation Methods:** Current Status & Modifications
- **Validation:** ^{56}Fe γ -ray Spectra Measurements
- **Conclusions & Future Work**



RPI Gaerttner Linear Accelerator (LINAC) Center

GAERTTNER LINAC LABORATORY
RENSSELAER POLYTECHNIC
INSTITUTE




Rensselaer

NAVAL NUCLEAR
LABORATORY



What do we measure?

Results of various neutron interactions



* indicates excited nucleus

$\nu_\gamma = \# \gamma$ -rays released per capture reaction

${}^A\text{X}$ and ${}^B\text{Y}$ are fission products

$\nu_n = \#$ neutrons released per fission reaction

Charged-particle and neutron-producing reactions can also occur but are not included in this research




Rensselaer

NAVAL NUCLEAR
LABORATORY

linac
The Gaertner LINAC Center

Thermal Neutron Capture Reactions



Project Goals

1

Develop **experimental methods** to measure **γ -ray cascades** produced by thermal neutron capture reactions

2

Update **simulation methods** to:

1. Accurately model **γ -ray cascades** emitted during compound nucleus de-excitation
2. Simulate **γ -ray cascades** travelling through detection systems

3

Use both methods to **assess the accuracy of thermal neutron capture induced γ -ray data** stored in nuclear data libraries

Why do we want to measure and simulate neutron capture γ -ray cascades?

Research Areas

- Nuclear structure studies
 - Energy levels, spins and parities
- Improving γ -ray de-excitation models

Applications

- γ -ray heating
- Radiation shielding
- Nuclear medicine
- Neutrino & particle physics
- Active neutron interrogation
 - Industry (coal, oil-well logging)
 - Controlled substance detection
 - Non-proliferation
 - Space exploration



Measurement Types

- Single Detector
- Coincidence Measurements
 - Multiple detectors to measure many γ rays emitted in a cascade
- **RPI capture γ -ray data can be used to test *both***

Gadolinium Neutron Capture Therapy (GdNCT)

- = tumor cell w/ neutron capturing isotope
- = Auger electron
- = IC electron

*figure not to scale



Active Neutron Interrogation

Experimental Methods


RPI Capture γ -Ray Multiplicity Detector

- **16 segment NaI(Tl) capture γ -ray multiplicity detector**

- 20 L of NaI(Tl) surrounding the sample
- A 1 cm thick B_4C ceramic sleeve (enriched to 99.5 atom% in ^{10}B) is used inside the detector to absorb neutrons scattering from the sample
- Up to 96% efficiency for detecting γ -ray cascades
- Located 25 m from the neutron-producing target
- Time-of-flight (TOF) method used to determine incident neutron energy

- **Used for neutron capture yield and γ -ray spectra measurements**

- Incident neutron energies: 0.001 eV – 3 keV



Simulation Methods

Step 1: Model Neutron Capture γ -Ray Cascades using DICEBOX

DICEBOX/ENSDF (+Firestone)


Models full γ -ray cascades using evaluated nuclear data (ENSDF + RIPL-3)

Input tuned to R. B. Firestone et. al., Phys. Rev. C **95**, 014328 (2017)

EGAF

Shows experimentally measured primary γ -ray lines (does not necessarily represent the full cascade)

$^{56}\text{Fe}(n,\gamma)$, $E_n = \text{thermal}$



Simulation Methods

Step 2: Transport γ -ray cascades through the RPI Capture Detector System

MCNP-6.2/ACE



- Extracts γ -ray data from ACE files (ENDF/B-VIII.0)
- Total energy deposition spectra is expected to disagree with experimental data because the simulation does not include coincidence events

MCNP-6.2/CGM

- Cascading γ -Ray Multiplicity
- Produces correlated secondary γ -ray emissions (cascades)
- γ rays transported through the detector geometry using MCNP-6.2

mod-MCNP-6.2/DICEBOX


- γ -ray cascades generated using DICEBOX
- γ -ray cascades transported through the detector geometry using a modified version of MCNP-6.2
- Generates an output file to tally γ -ray energy deposition in detector segments (enables event-by-event analysis including coincidence)



Validation of Experimental + Simulation Methods


Compare Experimental γ -ray Spectra to Simulated Results

$^{56}\text{Fe}(n,\gamma)$ Single Detector Spectra



*Black data points = average single detector response
Gray band indicates the range of the 16 detectors*


$^{56}\text{Fe}(n,\gamma)$ Total Energy Deposition (all 16 detectors) Spectra



Gray band indicates the uncertainty in total energy deposition (size is on the order of the data points)

Further Validation

^{59}Co Thermal Neutron Capture



DICEBOX/ENSDF


Models full γ -ray cascades using evaluated nuclear data (ENSDF)

EGAF


Experimentally measured γ -ray intensities

$^{59}\text{Co}(n,\gamma)$ Single Detector

Black data points = average single detector response
Gray band indicates the range of the 16 detectors



$^{59}\text{Co}(n,\gamma)$ Total Energy Deposition (all 16 detectors)




Rensselaer

NAVAL NUCLEAR
LABORATORY

linac
The Gaertner LINAC Center

Further Validation

^{55}Mn Thermal Neutron Capture



DICEBOX/ENSDF

Models full γ -ray cascades using evaluated nuclear data (ENSDF)

DICEBOX/ENSDF+levels


Models full γ -ray cascades using evaluated nuclear data (ENSDF)
+ additional levels that DICEBOX previously excluded

EGAF


Experimentally measured γ -ray intensities

$^{55}\text{Mn}(n,\gamma)$ Single Detector

*Black data points = average single detector response
Gray band indicates the range of the 16 detectors*



$^{55}\text{Mn}(n,\gamma)$ Total Energy Deposition (all 16 detectors)



Conclusions

- The experimental and simulation methods have been developed to test the accuracy of thermal neutron capture γ -ray data in nuclear data libraries.
- Measuring capture γ -ray spectra with the RPI Capture γ -Ray Multiplicity Detector system has been validated using the $^{56}\text{Fe}(n,\gamma)$ measurement.
 - **When the neutron capture γ -ray cascade data is well-known, the experimental γ -ray energy spectra can be accurately simulated using mod-MCNP-6.2/DICEBOX.**
- Validation of the system has been extended with ^{55}Mn and ^{59}Co thermal neutron capture measurements

Future Work

Complete the analysis of experimental capture γ -ray spectra for $^{\text{nat},235}\text{U}$ and compare to **mod-MCNP-6.2/DICEBOX** simulations

- Challenge: separate fission-induced γ rays from capture





Thank you! Questions?

This research was performed under appointment to the Rickover Fellowship Program in Nuclear Engineering sponsored by Naval Reactors (NR) Division of the National Nuclear Security Administration (NNSA).